




Transforming Patient Experience



SMART AND SIMPLE MEDICAL DEVICES PLATFORM FOR PREPARATION AND INJECTION




To learn more, visit www.eveon.eu
and contact us contact@eveon.eu

How to deliver complex drugs easily, automatically and safely?

New formulations, such as lyo products, highly-viscous products, suspensions, emulsions, radiopharmaceuticals, cell therapy, multidoses drugs, vaccines..., are challenging to deliver both for patients and healthcare professionals. EVEON offers modular solutions for smart and simple preparation and delivery.




Biologics



Highly viscous



Lyo
Multi-therapies



Suspensions
Emulsions



Time consuming
for the patient and
healthcare professionals



Imprecise
not accurate dosage



Complex formulations
lyo, viscous, suspensions,
emulsions, multitherapies,
radiopharmaceuticals,...



**Problematic of preparation
& delivery**
heterogeneous mixing, dosing
errors, user dependent



**Our custom-made solutions offer an
easy-to-use and patient compliant
medical device for complex therapies.**

Intuity® Ject, custom-made platform for smart and simple preparation



- Reusable device
- Disposable cassette adapted to all standard primary containers
1 cassette = 1 treatment
Adapted to mono & multi vials for SC or injection
- Connected & smart e-device
Data tracking, data sharing via GSM, drug and skin detection, automatized needle actuation

Specifications - Intuity® Ject

Micropumps technology → Motor driven oscillo rotary pump	
Displacement volumes	Custom and/or adapted for 300µL, 80µL, 20µL
Dosage	From 20µL (accuracy 7%) to several mL (accuracy 1%), up to 1L
Flowrate	0.4 µL/min to 90 mL/min with accuracy <2%
Pressure range	From -0.8bar to >10bar

Fluidic distribution and switching → 2 way valves (open/close) or multi-way distributors	
Tightness	5bar
Multiple positions pad/plan distributor	For complex fluids distribution to multiple containers
Ease of use	Minimum breakaway torque, actuating torque or linear force