



INNOVATION

Strong patent portfolio and know-how to secure your projects

OPERATIONAL EXCELLENCE
Part of JBT Hub Up alliance
securing your project from
conception to manufacturing



EXPERTISE

Cross-functional teams in fluidics, mechanics/ plasturgy, electronics/ software and QA



Our expertise enables us to adapt our QA/Regulatory system to your product regulatory status



AGILITY

Company size and project management system allow us to be agile and innovative problem solvers with a track record of providing outside-thebox solutions

SMART PREPARATION AND DELIVERY

TECHNOLOGY FOR BEAUTY



To learn more, visit www.eveon.eu and contact us contact@eveon.eu

How to do personalized preparation and delivery?

How to do personalized preparation and delivery?

The beauty market is focusing on customers' demand for personalized products and solutions. The technologies developed by EVEON meet this challenge and offer beauty tech devices with modular solutions for automated preparation and delivery.

Intuity® Mix, an extemporaneous & personalized mixing platform to prepare and deliver precise doses



Plug-and-play

Easy to use with only one user action From mixing to dispensing



Digital

Web app for users Personalization & Monitoring



Mixing solution

For a large range of preparations: emulsions, suspensions, gel/cream mixing...



Customization

Personalized formulation preparation on-demand

EVEON Intuity® Mix

- Reusable device with single-use disposable cassettes
- Adapted to any primary containers
- Mixing liquid/liquid (ie dilution, emulsification...) or liquid/solid (ie suspension, dispersion...) or gel mixing

Intuity® Spray, airless nozzle for a precise and controlled delivery of mist



Airless, propellant free

Sterility, eco-friendly

Fine and soft mist

Droplet diameter ~20µm

Low impact pressure for a

better user experience





Precise and personalized dosing

Controlled flow rate User-independent



Digital

Human-machine interface Personalization & Monitoring

EVEON spray nozzle

- For continuous delivery
- Accurate dosage
- From low to highly viscous formulations (to 8.000 Pa.s)
- Adapted to primary containers

Examples of targeted areas







Aesthetic F

Hair care

Dermatology